


CanDo® white foam rollers

high quality extruded foam roller

- extruded polyethylene 2.0 lb. density foam rollers provide moderate firmness
- shorter rollers save space in clinics and homes
- 6" diameter roller is the most popular
- use for positioning, balance, postural and muscle re-education, spinal stabilization, body awareness and coordination and strengthening activities
- use half-round roller for ankle stretching and knee rehabilitation; or as rocker board by standing on flat side
- foam rollers are strong and durable enough to be used by large adults and active children
- available in 12" to 48" lengths with either a round (circular) or half-round (semi-circular) profile
- can be used alone or with other exercisers such as CanDo® exercise band, tubing, weights and balls
- purchase by the case and save (in red below)



HALF-ROUND		ROUND	
8"	4"H x (8" Dia.)	8"	8" diameter
6"	3"H x (6" Dia.)	6"	6" diameter
4"	2"H x (4" Dia.)	4"	4" diameter



Jumbo™ 8" round

8" Jumbo™ white		
dia x L		
30-2260	8" x 36"	40.00
30-2261	8" x 12"	20.00



6" round

6" white		
dia x L		
* 30-2100	6" x 36"	20.00
* 30-2101	6" x 12"	10.00
**30-2100-12	6" x 36"	225.00
**30-2101-36	6" x 12"	342.50
30-2108	6" x 48"	28.00
30-2140	6" x 30"	18.00
30-2141	6" x 24"	16.00
30-2142	6" x 18"	14.00



4" round

4" white		
dia x L		
30-2102	4" x 36"	15.00
30-2103	4" x 12"	7.50



Slim™ 3" round

3" Slim™ white		
dia x L		
30-2106	3" x 36"	11.25
30-2107	3" x 12"	5.65



Jumbo™ 8" half-round

8" Jumbo™ white		
dia x L		
30-2270	8" x 36"	28.00
30-2271	8" x 12"	14.00



6" half-round

6" white		
dia x L		
* 30-2110	6" x 36"	14.00
* 30-2111	6" x 12"	7.00
**30-2110-24	6" x 36"	317.50
**30-2111-72	6" x 12"	477.50
30-2118	6" x 48"	19.60
30-2143	6" x 30"	12.60
30-2144	6" x 24"	11.20
30-2145	6" x 18"	9.80



4" half-round

4" white		
dia x L		
30-2112	4" x 36"	10.50
30-2113	4" x 12"	5.25



Slim™ 3" half-round

3" Slim™ white		
dia x L		
30-2116	3" x 36"	7.90
30-2117	3" x 12"	3.95

other sizes available on request

* most popular sizes

**buy by the case and save